

MPLC CHIRAL COLUMNS

# CHIRALFLASH<sup>®</sup> IA/IC/ID/IE/IF

Chiral TLC for method  
development of CHIRALFLASH

# 2L-ChiralTLC<sup>®</sup> IA/IC/ID/IE/IF



The new style of Chiral Separation



## MPLC CHIRAL COLUMNS

### High column efficiency

The 20 $\mu$ m CSP packed in CHIRALFLASH is specially designed. All solvents which are compatible with silica-gel based columns can be used.

### Easy to connect

To various MPLCs  
Coupling parts for all MPLCs are available

### Easy Method Development

From analytical HPLC to MPLC, the condition can be easily developed based on that of analytical HPLC column.

## 2L-ChiralTLC<sup>®</sup>

Achieved the UV detection by applying bilayer form which consisted of chiral layer and silica layer

By applying 2L-ChiralTLC, develop the preparative methods easily and quickly


We prepare IA/IC/ID/IE/IF types of 2L-ChiralTLC to correspond the series of CHIRALFLASH which commercially available

# CHIRALFLASH® IA/IC/ID/IE/IF

## Specifications


### CHIRALFLASH® IA

End cap : 1/4-28 UNF -Female  
 Chiral recognition : Amylose tris  
 reagent (3,5-dimethylphenylcarbamate)  
 Particle size : IA 20 μm




### CHIRALFLASH® IC

End cap : 1/4-28 UNF -Female  
 Chiral recognition : Cellulose tris  
 reagent (3,5-dichlorophenylcarbamate)  
 Particle size : IC 20 μm




### CHIRALFLASH® ID

End cap : 1/4-28 UNF -Female  
 Chiral recognition : Amylose tris  
 reagent (3-chlorophenylcarbamate)  
 Particle size : ID 20 μm




### CHIRALFLASH® IE

End cap : 1/4-28 UNF -Female  
 Chiral recognition : Amylose tris  
 reagent (3,5-dichlorophenylcarbamate)  
 Particle size : IE 20 μm



### CHIRALFLASH® IF

End cap : 1/4-28 UNF -Female  
 Chiral recognition : Amylose tris  
 reagent (3-chloro-4-methylphenylcarbamate)  
 Particle size : IF 20 μm




Column size	Packing size (mm) Tube size (mm)	30×100 38×150
Column materials	—	Fluoroplastic
Packing materials	CSP CSP amount (g)	IA 20μm, IC 20μm, ID 20μm, IE 20μm, IF 20μm 40
Bed volume	mL	50
Pressure limitation	MPa	Maximum pressure 1.5
Typical flow rate Sample load	mL/min.	12
Sample load	—	About 50mg ~ 100mg
Mobile phase	Alkanes, Alcohols, Ethyl acetate, THF, dichloromethane, etc.	
Column storage	If not using a column for one week or more, store the column in ethanol	



\* Important reminder

- Don't open the end-cap, otherwise it could negatively affect on separation.
- This column must be used for only investigational purpose.
- The maximum flow rate depends on the mobile phase viscosity (mobile phase composition), and should be adjusted in accordance with the upper limit of the column backpressure (i.e. 1.5 MPa).
- Prior to any preparative work, it is highly recommended that the CHIRALFLASH column be flushed with at least 300ml ~ 600 ml of eluent.
- The column may be cleaned by pumping solvent in the reverse direction (backflush) at half the normal flow rate for that solvent.
- Avoid strong mechanical shock to the column (for example dropping).

## Application data



# 2L-ChiralTLC® IA/IC/ID/IE/IF

## Specifications

Size (wide×length)	200mm×100mm
Layer Thickness	approximately 270μm
Particle size	CSP : 20μm
Recommended sample amount	1~5μL
Detection	Silica-gel layer includes a fluorescence indicator (UV254nm). Samples which have ultraviolet adsorption can be observed as spots (shadow).

※Important reminder

- The layer may be fallen off the plate depending on organic solvents.
- An impact shock or stress should be avoided as, the base material is aluminum. Especially when cutting the plate
- Don't press and rub the surface of TLC plates. Silica-gel layer may be fallen off. When you cut a TLC plate, a protection of the surface of the TLC plate is recommended.

### Chiral recognition reagent

### Sample spotting

- Spot the sample on "sample apply zone". ("Sample apply zone" is the bottom end of the TLC plate where only the CSP is coated) Recommended spotting position is within 10 mm of the bottom end of the TLC.
- Before spotting, it is recommended to verify the UV detection of the sample by spotting on the top of the silica-gel layer zone.

### If sample spots can't be observed

Spray and dry an organic solvent (e.g. Ethanol) on the surface of TLC. It may help the detection  
An accompanying atomizer is for this propose.

### For acidic samples or basic samples

For basic samples or acidic samples, it is necessary to add an additive into the developing solvent in order to get appropriate spots, otherwise broad and/or tailing spot form may be observed.

Please follow the bellow procedure.






1. Add 0.1 volume percent of acid (for acidic samples) or base (for basic samples) to developing solvent.
2. Before developing, soak TLC plate for about 1 minute in developing solvent and dry it.
3. Operate developing with the TLC plate prepared in "2" and solvent of "1".

### For the samples which don't have enough UV adsorption

Use coloring reagent to detect the spot. (Coloring reagent : phosphomolybdic acid-ethanol, iodine stain, p-anisaldehyde, or ninhydrin)

If you need to use coloring reagent, it is recommended to verify the coloring and the detection of the sample by spotting the sample on a commercially available silica-based TLC plate before using 2L-ChiralTLC


## Application data

Product name	2L-ChiralTLC® IA	2L-ChiralTLC® IC	2L-ChiralTLC® ID	2L-ChiralTLC® IE	2L-ChiralTLC® IF
Sample name	Warfarin Sodium	4-Benzoyloxy-2-azetidinone	Aminogluthethimide	Wieland Miescher ketone (WMK)	Furoin
Solvent	n-Hex/2-Propanol/TFA =50/50/0.1(v/v/v)	n-Hex/2-Propanol =50/50(v/v)	AE/DEA=100/0.1(v/v)	n-Hex/EtOH=60/40 (v/v) (Two times development)	n-Hex/EtOH/MeOH =75/15/15(v/v/v)
Sample conc.	25,000ppm	100,000ppm	20,000ppm	10,000ppm	50,000ppm
Sample amount	1μL	1μL	1μL	1μL	1μL
Detection	UV 254nm	UV 254nm	UV 254nm	UV 254nm	UV 254nm
	 Rf <sub>1</sub> =0.78 Rf <sub>2</sub> =0.60	 Rf <sub>1</sub> =0.41 Rf <sub>2</sub> =0.28	 Rf <sub>1</sub> =0.44 Rf <sub>2</sub> =0.16	 Rf <sub>1</sub> =0.51 Rf <sub>2</sub> =0.41	 Rf <sub>1</sub> =0.28 Rf <sub>2</sub> =0.16

※TFA=trifluoroacetic acid, DEA=diethylamine, AE=ethyl acetate


## Method development to MPLC

### Flavanone




TLC : 2L-ChiralTLC® IA  
Solvent : n-Hex/EtOH=90/10(v/v)  
Sample conc. : 5,000ppm  
Sample amount : 1 μL  
Detection : UV 254nm

Mobile phase : CHIRALFLASH® IA  
Flow rate : n-Hex/EtOH=90/10(v/v)  
Sample conc. : 12mL/min.  
Sample conc. : 30g/L  
Injection volume : 30, 60, 90mg  
Detection : UV 254nm




## MPLC CHIRAL COLUMNS

Column name	I.D.(mm)	Length(mm)	Particle size (μm)	CSP amount (g)	Code
CHIRALFLASH® IA	30	100	20	40	80M73
CHIRALFLASH® IC	30	100	20	40	83M73
CHIRALFLASH® ID	30	100	20	40	84M73
CHIRALFLASH® IE	30	100	20	40	85M73
CHIRALFLASH® IF	30	100	20	40	86M73



## Analytical HPLC columns for MPLC method development

Column name	I.D.(mm)	Length(mm)	Particle size (μm)	Code
CHIRALPAK® IA (20 μm)	4.6	100	20	80223
CHIRALPAK® IC (20 μm)	4.6	100	20	83223
CHIRALPAK® ID (20 μm)	4.6	100	20	84223
CHIRALPAK® IE (20 μm)	4.6	100	20	85223
CHIRALPAK® IF (20 μm)	4.6	100	20	86223



## Injection Columns

Column name	I.D.(mm)	Length(mm)	Number of columns	Code
Injection column S	15	44	10	00M01
Injection column M	20	75	10	00M02
Injection column L	26	80	10	00M03

## Column Adaptor For Injection Columns


Items	Code
Column Adaptor S	00M10
Column Adaptor M	00M11
Column Adaptor L	00M12

## Column joint for CHIRALFLASH®

Items	Specification	Code
Column Rack Adapter	1/4-28UNF	00M04
Injection Adapter	1/4-28UNF	00M05
Lure Lock Adapter(Male)	1/4-28UNF	00M06
Syringe Union	1/4-28UNF	00M07
"LL Adapter [INTRM Connector -Female]"	1/4-28UNF	00M08
"LL Adapter [INTRM Connector -Male]"	1/4-28UNF	00M09


## Column Holder For Injection Columns

Items	Code
Injection columns S	00M13
Injection columns M	00M14
Injection columns L	00M15



## 2L-ChiralTLC®

Items	Sheet size (mm)	Number of sheets	Code
2L-ChiralTLC® 5 types	200×100	5	5T5A
2L-ChiralTLC® IA	200×100	2	80T2A
2L-ChiralTLC® IC	200×100	2	83T2A
2L-ChiralTLC® ID	200×100	2	84T2A
2L-ChiralTLC® IE	200×100	2	85T2A
2L-ChiralTLC® IF	200×100	2	86T2A



※ 5 types includes IA, IC, ID, IE, IF.

## DAICEL CORPORATION [TOKYO HEAD OFFICE]

JR Shinagawa East Bldg. 2-18-1, Konan, Minato-ku, Tokyo 180-8230, Japan  
 Phone: +81-3-6711-8222(Direct)  
 Facsimile: +81-3-6711-8228  
 URL: <http://www.daicelchiral.com/>  
 E-mail: [chiral@jp.daicel.com](mailto:chiral@jp.daicel.com)

## DAICEL CHIRAL TECHNOLOGIES (INDIA) PVT. LTD

Lab No.4A, Phase III, ICICI Knowledge Park Genome Valley, Turkpally,  
 Shameerpet, Range Reddy Dist, Hyderabad-500 078, A.P., Republic of India  
 Phone: +91-40-23480103  
 Facsimile: +91-40-23480104